

# Blob\_upload

April 27, 2016

## 1 MDCS REST api demo – data from Diffusion couple

```
In [1]: import pip
```

```
requires = ['numpy', 'xmldict']

installed_packages = pip.get_installed_distributions()
installed_packages_list = sorted(["%s==%s" % (i.key, i.version) for i in installed_packages])

matching = [[libs for libs in installed_packages_list if x in libs] for x in requires]

if len(matching) == len(requires):
    print "requirement fullfilled"
else:
    print "something is missing"
```

requirement fullfilled

NOTE: if the check is failed, execute the following cell

```
In [ ]: #!/conda install numpy
        #!/pip install xmldict
```

```
In [2]: import os
```

```
file_check = ['DataModelDict.py', 'GE-DiffusionCouple-IN100-IN718.xml', 'GE-DiffusionCouple-IN718',
              'GE-DiffusionCouple-R88-R95.xml', 'mdcs', 'data', 'fig']

files = os.listdir(".")
matching = filter(None, [[fs for fs in files if x in fs] for x in file_check])

if len(matching) == len(file_check):
    print "files in the folder"
else:
    print "something is missing"
```

files in the folder

```
In [3]: # standard python libraries
import numpy as np
import glob
import requests
import xmldict
```


```
# ----- Matplotlib for plotting
# NEED to be installed (anaconda)

# http://matplotlib.org/
import matplotlib.pyplot as plt
import matplotlib.colorbar as clb
import matplotlib.patches as mpatches
import matplotlib.lines as mlines
import matplotlib.gridspec as gridspec
from matplotlib import ticker
from matplotlib.ticker import LogLocator
%matplotlib inline

# customized library by Zachary Trautt
import mdcs
# customized library by Lucas Hale
from DataModelDict import DataModelDict
```

```
In [4]: from IPython.display import Image
Image("fig/file_problem.png", width=600)
```

Out[4]:



## 1.1 How to manage .xml and image files? – the update files are prepared in the folder “data”

### 1.1.1 steps:

#### 1.1.2 (1) upload the draft XMLs

#### 1.1.3 (2) upload image files

#### 1.1.4 (3) query the data in MDCS

#### 1.1.5 (4) update the content with image URL

#### 1.1.6 (5) delete the draft XMLs from MDCS

#### 1.1.7 (6) upload the updated XMLs

```
In [5]: # ---- setup user information to sign in MDCS
        USER = "admin"
        PSWD = "admin"
        MDCS_URL = "http://127.0.0.1:8000"
```

### 1.1.8 obtain the file names

```
In [7]: xml_files = glob.glob("data/*.xml")
        print "Files: ", len(xml_files), xml_files

        filetitle = []
        xmlfile = []
        fig_files = []

        for i in range(len(xml_files)):
            filetitle.append(xml_files[i][5:-4])
            fig_files.append(str(xml_files[i][0:len(xml_files[i])-3] + '.tif'))

        print "XML files: ", filetitle
        print "Image files: ", fig_files
```

```
Files: 3 ['data/GE-DiffusionCouple-R88-R95.xml', 'data/GE-DiffusionCouple-IN718-R95.xml', 'data/GE-Di-
XML files: ['GE-DiffusionCouple-R88-R95', 'GE-DiffusionCouple-IN718-R95', 'GE-DiffusionCouple-IN100-IN
Image files: ['data/GE-DiffusionCouple-R88-R95.tif', 'data/GE-DiffusionCouple-IN718-R95.tif', 'data/GE-
```

### 1.1.9 Using blob function to upload the image files and return the URL

```
In [8]: figURL = []

        for i in range(len(xml_files)):
            figURL.append(mdc.blobs.upload(fig_files[i], MDCS_URL, USER, PSWD, 'cert'))

        print figURL

[u'http://127.0.0.1:8000/rest/blob?id=5720f57a1ff0f31620f6cc42', u'http://127.0.0.1:8000/rest/blob?id=5'
```

### 1.1.10 upload the XML files

make sure the template\_title is correct on your MDCS instance

```
In [9]: for i in range(len(xml_files)):
        print mdc.curate_as(xml_files[i], filetitle[i], MDCS_URL, \
            USER, PSWD, 'cert', template_title="mod_demo_diffusion0224")
```

```
created
created
created
```

1.1.11 demo the function “query” to search the data by filetitles in database:

1.1.12 (1) the data content is stored in python list “qres”

1.1.13 (2) the specific “id” assigned by MDCS is kept in file\_id

```
In [10]: qres = []
         file_id = []

         for i in range(len(xml_files)):
             content = mdcs.explore.query(MDCS_URL,USER,PSWD,'cert', \
                                           None,str({'title':filetitle[i]}))
             qres.append(xmltodict.parse(content[0]['content']))
             file_id.append(content[0]['_id'])
```

1.1.14 update XML files with image URL

```
In [11]: res = []

         for i in range(len(xml_files)):
             qres[i]['experiment']['experimentType']['chemicalDiffusivity'] \
             ['dataAnalysis']['micrographOfStructure']['fileExtension'] = figURL[i]
             res.append(xmltodict.unparse(qres[i]))
```

## 1.2 output the content as new XMLs

1.2.1 WARNING: Occassionally, Jupyter may cause issues to output large file. You may try to run this cell multiple times.

```
In [12]: filename = []

         # write information to interdiffusion.xml
         for i in range(len(xml_files)):
             filename.append(str(filetitle[i] + str('.xml')))
             diffusion = open(filename[i], 'w')
             diffusion.write(res[0].encode('utf8'))
             diffusion.close()
```

## 1.3 delete the old XMLs in MDCS

```
In [13]: for i in range(len(xml_files)):
         mdcs.explore.delete(file_id[i],MDCS_URL,USER,PSWD,cert=None)
```

## 1.4 upload the revised XMLs

```
In [14]: for i in range(len(xml_files)):
         print filename[i],mdcs.curate_as(filename[i],filetitle[i]
         ,MDCS_URL,USER,PSWD,'cert',template_title="mod_demo_diffusion0224")
```

```
GE-DiffusionCouple-R88-R95.xml created
GE-DiffusionCouple-IN718-R95.xml created
GE-DiffusionCouple-IN100-IN718.xml created
```

## 1.5 Visualize data

### 1.5.1 query, load, and plot the data

```
In [15]: content = mdcs.explore.query(MDCS_URL,USER,PSWD,'cert',
                                     None,str({'title':filetitle[0]}))
        qdata = DataModelDict(content[0]['content'])

        table = qdata.find('table')

In [16]: distable = []

        for row in table['rows'].iterlist('row'):
            disrow = []
            for column in row.iterlist('column'):
                disrow.append(column['#text'])

            distable.append(disrow)


        del distable[0]
```

### 1.6 Plot the data

```
In [17]: distable = np.array(distable)

        fig, ax = plt.subplots(figsize=(9, 6))
        for i in range(2,len(distable[0])):
            ax.plot(distable[:,0],distable[:,i],'o-')

        fig.tight_layout()
        plt.show()
```



```
In [ ]:
```